




ROLLOVER ALARM™

Instructions



- Plug power cord and sensor into the RollOver Alarm™. Plug the alarm into a power source.
- Place jersey headband around head. Tuck the sensor under the band at the back of the head. When lying in the face down position, the sensor should be parallel to the bed with the white dot facing the ceiling.
- Press the On/Off button. The alarm will sound until the Reset button is pressed. This is to indicate the unit is operational.
- If the head tilts or turns to a 30 degree angle for more than 8 seconds, the alarm will sound.
- To reset the alarm, return head to proper face down positioning and press the Reset button. The alarm will NOT reset until the head has been returned to the proper position.